

TL 1000 & TL 1000 ASOPS REPETITION RATE STABILIZATION WELL-TIMED FEMTOSECOND PULSES

Novanta develops high-precision lasers and essential accessories for cutting-edge ultrafast scientific and academic research through our globally recognized brands like Laser Quantum. Designed and built for accurate results, reliable operation, and straight-forward integration, our ultrafast lasers can be found in research and development facilities around the globe.

REPETITION RATE STABILIZATION


Novanta offers the timing stabilisation units TL-1000 and TL-1000-ASOPS as accessories to the its series of femtosecond lasers. The TL-1000 allows the tight phase-lock of an oscillator repetition rate to an external reference such as a synthesiser or another mode locked laser. TL-1000-ASOPS enables a repetition rate offset-lock between two femtosecond oscillators (e.g. two taccor lasers or two gecco lasers) at a repetition rate difference between 20 Hz and 20 kHz. It also permits synchronisation of two lasers with a jitter below 60 fs (synch. mode). System parameters are accessible via touchscreen and USB port

Repetition rate stabilisation is essential for applications requiring a well-defined timing relation between a femtosecond oscillator and a reference signal. Such applications are for example two-color pump-probe spectroscopy using two synchronised mode locked lasers, optical experiments synchronised to a pulsed electron source, or high speed asynchronous optical sampling (ASOPS) experiments. ASOPS is an ultra-rapid and precise time-domain spectroscopy technique pioneered by Laser Quantum, outperforming classical approaches by orders of magnitude in measurement speed and noise performance.

The TL-1000 and TL-1000-ASOPS are specifically designed to serve these applications. Repetition rate fluctuations of a free-running oscillator are efficiently suppressed. TL-1000-ASOPS permits high-speed ASOPS based ultrafast time-domain spectroscopy using two femtosecond oscillators with >60 fs time-delay resolution.


OPTIONAL FEATURE & FUNCTIONALITY

Optional ASOPS experiment support - the TL-1000-ASOPS is capable of driving and monitoring a photoreceiver and a Tera-SED THz emitter element via its umbilical port. This feature is beneficial for high-speed ASOPS experiments.



Functionality


- Repetition rate stabilisation
- Repetition rate offset locking
- Low residual timing jitter
- ASynchronous Optical Sampling (ASOPS)
- Time Domain Spectroscopy



The TL-1000-ASOPS controls the time delay between pump and probe pulses by slightly differing repetition rates ($f_R \sim 1$ GHz). The difference in Δf_R determines the scan rate, while the measurement window is given by the inverse of the repetition rate $1/f_R = 1$ ns. Here, Δf_R is 10 kHz, i.e. the time-delay is repetitively ramped with a 100 μ s period.

TL-1000 & TL-1000 ASOPS SPECIFICATIONS

DIMENSIONS (mm)



SPECIFICATIONS

Specification*	TL-1000	TL-1000-ASOPS
Repetition Rate range ¹	75 MHz to 10 GHz	79.68 MHz to 80.32 MHz 0.996 GHz to 1.004 GHz
Timing jitter ^{2,3} : taccor models	≤ 100 fs (0.1 Hz to 100 kHz)	n/a
Timing jitter ^{2,3} : venteon/gecco models	≤ 300 fs (0.1 Hz to 100 kHz)	n/a
Timing jitter in sync mode ³	n/a	< 60 fs (0.1 Hz to 100 kHz) @ 1 GHz < 100 fs (0.1 Hz to 100 kHz) @ 80 MHz
Repetition Rate offset	n/a	0 (synch. mode) or 20 Hz to 20 kHz
Time resolution in ASOPS	n/a	< 60 fs @ 1 GHz < 100 fs @ 80 MHz
Trigger Signal	n/a	TTL level at offset frequency, ≤ 10 ns rise time
Power Requirements	110 or 220 VAC (60/50 Hz)	110 or 220 VAC (60/50 Hz)

Notes:

1 - Ranges must be selected upon order.

2 - Relative to customer provided signal at +7 dBm (50 Ohm) with <-125 dBc/Hz phase-noise above 10 kHz offset from carrier.

3 - If used with a Laser Quantum femtosecond oscillator and suitable built-in piezo supported cavity mirror(s).

CONTACT US

Americas, Asia Pacific

Novanta Headquarters
Bedford, USA
P +1-781-266-5700

Photonics@Novanta.com

Europe, Middle East, Africa

Stockport, United Kingdom
P +44-161-975-5300

Novanta Europe GmbH
Wackersdorf, Germany
P +49 9431 7984-0

Milan, Italy
P +39-039-793-710

Photonics@Novanta.com

China

Novanta Sales & Service Office
Shenzhen, China
P +86-755-8280-5395

Suzhou, China
P +86-512-6283-7080

Photonics.China@Novanta.com

Japan

Novanta Service & Sales Office
Tokyo, Japan
P +81-3-5753-2460

Photonics.Japan@Novanta.com